

Product datasheet for **TP325217L**

CSRP3 (NM_001127656) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cysteine and glycine-rich protein 3 (cardiac LIM protein) (CSRP3), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225217 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPNWGGGAKCGACEKTVYHAEIQCNGRSFHKTCFHCMACRKALDSTTVAAHESIYCKVCYGRRYGPKG IGYGQGAGCLSTDTGEHLGLQFQQSPKPARSVTTSNPSKFTAKFGESEKPCRCGKSVYAAEKVMGGGKPW HKTCFRCAICGKSLESTNVTDKDGELYCKVCYAKNFGPTGIGFGGLTQQVEKKE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	20.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Enzyme substrate (PMID: 27281159) Binding assay (PMID: 27281159) Co-immunoprecipitation (PMID: 27281159)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001121128




[View online »](#)


Locus ID: 8048
UniProt ID: [P50461](#)
RefSeq Size: 1464
Cytogenetics: 11p15.1
RefSeq ORF: 582
Synonyms: CLP; CMD1M; CMH12; CRP3; LMO4; MLP

Summary: This gene encodes a member of the CSRP family of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this protein is found in a group of proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Mutations in this gene are thought to cause heritable forms of hypertrophic cardiomyopathy (HCM) and dilated cardiomyopathy (DCM) in humans. Alternatively spliced transcript variants with different 5' UTR, but encoding the same protein, have been found for this gene. [provided by RefSeq, Jul 2008]


Product images:




CSRP3 is the substrate of O-glycosidase, according to Western blot (WB) and *Amaranthus caudatus* lectin blot (LB) analyses of recombinant human CSRP3. Recombinant proteins expressed in *E. coli* (unglycosylated negative control) and HEK293 cells (potentially glycosylated reference) (OriGene [TP325217]) were analyzed after treatment with sialidase and O-glycosidase. Figure cited from PLoS ONE, PMID: 27281159



CSRP3 binds to *Amaranthus caudatus* lectin (ACA), according to Western blot (WB) and ACA lectin blot (LB) analyses of recombinant human CSRP3. Recombinant proteins expressed in *E. coli* (unglycosylated negative control) and HEK293 cells (potentially glycosylated reference) (OriGene [TP325217]) were analyzed after treatment with sialidase and O-glycosidase. Figure cited from PLoS ONE, PMID: 27281159



Altered O-glycosylation on CSRP3 in the left ventricle of Dahl salt-sensitive rats, according to Western blot (WB) and *Amaranthus caudatus* lectin blot (LB) analyses of CSRP3 from the left ventricle extracts of Dahl salt-sensitive rats under a low-salt (LS) or high-salt (HS) diet for 12 or 16 weeks. CSRP3 in the extracts was immunoprecipitated, denatured, separated by SDS-PAGE, and analyzed. Recombinant human CSRP3 (OriGene [TP325217]) served as a control for immunoprecipitation with the anti-CSRP3 antibody (+) or normal IgG (-). Figure cited from PLoS ONE, PMID: 27281159



Coomassie blue staining of purified CSRP3 protein (Cat# [TP325217]). The protein was produced from HEK293T cells transfected with CSRP3 cDNA clone (Cat# [RC225217]) using MegaTran 2.0 (Cat# [TT210002]).