

Product datasheet for **TP325251L**

CLDN19 (NM_001123395) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human claudin 19 (CLDN19), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225251 representing NM_001123395 Red =Cloning site Green =Tags(s)
	<p>MANSGLQLLGYFLALGGWVGIIASTALPQWKQSSYAGDAITAVGLYEGLWMSCASQSTGQVQCKLYDSL LALDGHIQSARALMVAVLLGFVAMVLSVGMKCTRVGDSNPIAKGRVAIAGGALFILAGLCTLTAVSWY ATLVTQEFFNPSTPVNARYEFGPALFVGWASAGLAVLGGSFLLCCTCPEPERPNSSQPYPYRPGPSAAAREY V</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	21.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001116867</u>
Locus ID:	149461
UniProt ID:	<u>Q8N6F1</u>



[View online »](#)

Cytogenetics: 1p34.2

RefSeq ORF: 633


Synonyms: HOMG5

Summary: The product of this gene belongs to the claudin family. It plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. Defects in this gene are the cause of hypomagnesemia renal with ocular involvement (HOMGO). HOMGO is a progressive renal disease characterized by primary renal magnesium wasting with hypomagnesemia, hypercalciuria and nephrocalcinosis associated with severe ocular abnormalities such as bilateral chorioretinal scars, macular colobomata, significant myopia and nystagmus. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jun 2010]

Protein Families: Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction

Product images:



Coomassie blue staining of purified CLDN19 protein (Cat# [TP325251]). The protein was produced from HEK293T cells transfected with CLDN19 cDNA clone (Cat# [RC225251]) using MegaTran 2.0 (Cat# [TT210002]).