

Product datasheet for **TP325272L**

FAM164B (ZC2HC1B) (NM_001013623) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 164, member B (FAM164B), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225272 representing NM_001013623 Red =Cloning site Green =Tags(s)
	 MAGAEPFLADGNQELFPCEVCGRRFAADVLERHGPICKKLFNRKRKPFSSLKQRLQGTDIPTVKKTPQSK SPPVRKSNWRQQHEDFINAIRSAKQCMLAIKEGRPLPPPPPSLNPDIYIQRPYCMRRFNESAAERHTNFC KDQSSRRVFNPAQTAAKLASRAQGRAQMGPKEPTVTSAVGALLQNRVLVATNEVPTKSLAMD PASGAK LRQGFSSKSKD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001013645</u>
Locus ID:	153918



[View online »](#)

UniProt ID: [Q5TFG8](#)
Cytogenetics: 6q24.2
RefSeq ORF: 666
Synonyms: C6orf94; dj468K18.5; FAM164B

Product images:



Coomassie blue staining of purified ZC2HC1B protein (Cat# [TP325272]). The protein was produced from HEK293T cells transfected with ZC2HC1B cDNA clone (Cat# [RC225272]) using MegaTran 2.0 (Cat# [TT210002]).