

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP325345

Bcl2 Binding component 3 (BBC3) (NM_001127240) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BCL2 binding component 3 (BBC3), transcript variant 1, 20 μg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225345 representing NM_001127240 <mark>Red</mark> =Cloning site Green=Tags(s)
•	MKFGMGSAQACPCQVPRAASTTWVPCQICGPRERHGPRTPGGQLPGARRGPGPRRPAPLPARPPGALGSV LRPLRARPGCRPRRPHPAARCLPLRPHRPTRRHRRPGGFPLAWGSPQPAPRPAPGRSSALALAGGAAPGV ARAQRPGGSGGRSHPGGPGSPRGGGTVGPGDRGPAAADGGRPQRTVRAAETRGAAAAPPLTLEGPVQSHH GTPALTQGPQSPRDGAQLGACTRPVDVRDSGGRPLPPPDTLASAGDFLCTM
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	26.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP 001120712</u>
Locus ID:	27113
UniProt ID:	<u>Q9BXH1, Q96PG8</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Bcl2 Binding component 3 (BBC3) (NM_001127240) Human Recombinant Protein – TP325345
Cytogenetics:	19q13.32
RefSeq ORF:	783
Synonyms:	JFY-1; JFY1; PUMA
Summary:	This gene encodes a member of the BCL-2 family of proteins. This family member belongs to the BH3-only pro-apoptotic subclass. The protein cooperates with direct activator proteins to induce mitochondrial outer membrane permeabilization and apoptosis. It can bind to anti-apoptotic Bcl-2 family members to induce mitochondrial dysfunction and caspase activation. Because of its pro-apoptotic role, this gene is a potential drug target for cancer therapy and for tissue injury. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2011]
Protein Families	Druggable Genome
Protein Pathway	s: Huntington's disease, p53 signaling pathway

Product images:



Coomassie blue staining of purified BBC3 protein (Cat# TP325345). The protein was produced from HEK293T cells transfected with BBC3 cDNA clone (Cat# [RC225345]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US