

## Product datasheet for **TP325413**

### PRR32 (NM\_001122716) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome X open reading frame 64 (CXorf64), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225413 representing NM_001122716 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MACIENVLGGHAPSPLVWSVDKNGNQLHHDMPLQCLSSKPEDDAEPWGQPQVPLRPSVNVLTDLDSKQL EWPSERTGSCIPLHSLRAHRHPYGPPIPAVAEESLATAEVNSSDALAGWRQEGQDAINVSWEVSGGPPALI VGGTKVNNGGTERGSNNARLHVALPQKGFFPPRGPQVRGPSHIPTLRSGIVMEVPPGNTRIACRGKLAH VSFPLRGPCHPMHNWPRPIPLSSSTPGLPSCSTVHCIFPPRPPIFNPFILTMPLPFAPPPIFGPPLPSYFA HFHSGGMPAPASPNREHS</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	31.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001116188</a></u>
Locus ID:	100130613



[View online »](#)

UniProt ID: [B1ATL7](#)  
 Cytogenetics: Xq25  
 RefSeq ORF: 894  
 Synonyms: CXorf64

### Product images:



Coomassie blue staining of purified PRR32 protein (Cat# TP325413). The protein was produced from HEK293T cells transfected with PRR32 cDNA clone (Cat# [RC225413]) using MegaTran 2.0 (Cat# [TT210002]).