

## Product datasheet for **TP325416L**

### CEACAM19 (NM\_001127893) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human carcinoembryonic antigen-related cell adhesion molecule 19 (CEACAM19), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225416 representing NM_001127893 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MEIPMGTQGCFSKLLLSASILVLWMLQGSQAALYIQKIPEQPQKNQDLLLSVQGVPTDFQDFNWYLGEE  
TYGGTRLFTYIPGIQRPQRDGSAMGQRDIVGFNGSMMLRRAQPTDSGTQVAITINSEWTMKAKTEVQV  
AEKNKELPSTHLPTNAGILAATIIGSLAAGALLISCIAYLLVTRNWRGQSHRLPAPRGQGSLSILCSAVS  
PVPSVTPSTWMATTEKPELGAHDAGDNNIYEVMPSPVLLVSPISDTRSINPARPLTPPHLQAEPENHQ  
YQDLLNPDPAPYCQLVPTS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	32.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001121365</a>
Locus ID:	56971



[View online »](#)

UniProt ID: [Q7Z692](#)  
Cytogenetics: 19q13.31  
RefSeq ORF: 897  
Synonyms: CEACM19; CEAL1  
Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified CEACAM19 protein (Cat# [TP325416]). The protein was produced from HEK293T cells transfected with CEACAM19 cDNA clone (Cat# [RC225416]) using MegaTran 2.0 (Cat# [TT210002]).