

Product datasheet for **TP325490M**

MSF (SEPT9) (NM_001113496) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens septin 9 (SEPT9), transcript variant 7, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225490 representing NM_001113496 Red =Cloning site Green =Tags(s)
	<p>MADTPRDAGLKQAPASRNEKAPVDFGYVGIDSILEQMRRKAMKQGFEFNIMVVGQSLGKSTLINTLFKS KISRKSVQPTSEERIPKTIEIKSITHDIEEKGVRMCLTVIDTPGFGDHINNENCWQPIMKFINDQYEKYL QEEVNINRKKRIPDTRVHCCLYFIPATGHSLRPLDIEFMKRLSKVWNIVPIAKADTLLEERVHFKQRI TADLLSNGIDVYPQKEFDESEDRLVNEKFREMIPFAVVGSDHEYQVNGKRILGRKTKWGTIEVENTTHC EFAYLRDLLIRTHMQNIKIDITSSIHFEAYRVKRLNEGSSAMANGMEEKEPEPEM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	38.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001106968
Locus ID:	10801



[View online »](#)

UniProt ID: [Q9UHD8](#)

Cytogenetics: 17q25.3


RefSeq ORF: 1005

Synonyms: AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1

Summary: This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2009]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified SEPTIN9 protein (Cat# [TP325490]). The protein was produced from HEK293T cells transfected with SEPTIN9 cDNA clone (Cat# [RC225490]) using MegaTran 2.0 (Cat# [TT210002]).