

Product datasheet for TP325494

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

Activin A Receptor Type IC (ACVR1C) (NM 001111033) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens activin A receptor, type IC (ACVR1C), transcript

variant 4, 20 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC225494 representing NM 001111033

or AA Sequence: Red=Cloning site Green=Tags(s)

MTRALCSALRQALLLLAAAAELSPGLKCVCLLCDSSNFTCQTEGACWASVMLTNGKEQVIKSCVSLPELN AQVFCHSSNNVTKTECCFTDFCNNITLHLPTDNGTWTQLWLVSEYHEQGSLYDYLNRNIVTVAGMIKLAL SIASGLAHLHMEIVGTQGKPAIAHRDIKSKNILVKKCETCAIADLGLAVKHDSILNTIDIPQNPKVGTKR YMAPEMLDDTMNVNIFESFKRADIYSVGLVYWEIARRCSVGGIVEEYQLPYYDMVPSDPSIEEMRKVVCD

QKFRPSIPNQWQSCEALRVMGRIMRECWYANGAARLTALRIKKTISQLCVKEDCKA

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Tag: C-Myc/DDK

Predicted MW: 37.3 kDa

Concentration: $>0.05 \mu g/\mu L$ as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeg: NP 001104503

Locus ID: 130399



UniProt ID: Q8NER5

Cytogenetics: 2q24.1 RefSeq ORF: 1008

Synonyms: ACVRLK7; ALK7

Summary: ACVR1C is a type I receptor for the TGFB (see MIM 190180) family of signaling molecules.

Upon ligand binding, type I receptors phosphorylate cytoplasmic SMAD transcription factors, which then translocate to the nucleus and interact directly with DNA or in complex with other transcription factors (Bondestam et al., 2001 [PubMed 12063393]).[supplied by OMIM, Mar


2008]

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Adherens junction, Chronic myeloid leukemia, Colorectal cancer, Endocytosis, MAPK signaling

pathway, Pancreatic cancer, Pathways in cancer, TGF-beta signaling pathway

Product images:



Coomassie blue staining of purified ACVR1C protein (Cat# TP325494). The protein was produced from HEK293T cells transfected with ACVR1C cDNA clone (Cat# [RC225494]) using MegaTran 2.0 (Cat# [TT210002]).