

Product datasheet for **TP325494L**

Activin A Receptor Type IC (ACVR1C) (NM_001111033) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens activin A receptor, type IC (ACVR1C), transcript variant 4, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225494 representing NM_001111033 Red =Cloning site Green =Tags(s)

MTRALCSALRQALLLLAAAAELSPGLKCVCLLCDSSNFTCQTEGACWASVMLTNGKEQVIKSCVSLPELN
AQVFCHSSNNVTKTECCFTDFCENNITLHLPTDNGTWTQLWLWSEYHEQGSLYDYLNRNIVTVAGMIKLAL
SIASGLAHLHMEIVGTQGKPAIAHRDIKSKNILVKKCETCAIADLGLAVKHDSILNTIDIPQNPKVGTKR
YMAPEMLDDTMNVNIFESFKRADIYSVGLVYWEIARRCSVGGIVEEYQLPPYDMVPSDPSIEEMRKVVCD
QKFRPSIPNQWQSCEALRVMGRIMRECWYANGAARLTALRIKKTISQLCVKEDCKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	37.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001104503
Locus ID:	130399



[View online »](#)

UniProt ID: [Q8NER5](#)

Cytogenetics: 2q24.1

RefSeq ORF: 1008


Synonyms: ACVRLK7; ALK7

Summary: ACVR1C is a type I receptor for the TGFB (see MIM 190180) family of signaling molecules. Upon ligand binding, type I receptors phosphorylate cytoplasmic SMAD transcription factors, which then translocate to the nucleus and interact directly with DNA or in complex with other transcription factors (Bondestam et al., 2001 [PubMed 12063393]).[supplied by OMIM, Mar 2008]

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Adherens junction, Chronic myeloid leukemia, Colorectal cancer, Endocytosis, MAPK signaling pathway, Pancreatic cancer, Pathways in cancer, TGF-beta signaling pathway

Product images:



Coomassie blue staining of purified ACVR1C protein (Cat# [TP325494]). The protein was produced from HEK293T cells transfected with ACVR1C cDNA clone (Cat# [RC225494]) using MegaTran 2.0 (Cat# [TT210002]).