

Product datasheet for TP325497

FBP1 (NM_001127628) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fructose-1,6-bisphosphatase 1 (FBP1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225497 protein sequence Red =Cloning site Green =Tags(s)

MADQAPFDTDVNTLTRFVMEGRKARGTGELTQLLSLCTAVKAISSAVRKAGIAHLYGIAGSTNVTGDQ
 VKKLDVLSNDLVMNMLKSSFATCVLVSEEDKHAIIVEPEKRGKYVVCFDPLDGSSNIDCLVSVGTIFGIY
 RKKSTDEPSEKDALQPGRNLVAAGYALYGSATMLVLAMDCGVNCFMLDPAIGEFILVDKDVKIKKKGKIY
 SLNEGYARDFDPAVTEYIQRKKFPPDNSAPYGARYVGSVMADVHRTLTVGGIFLYPANKKSPNGKLRLLY
 ECNPMAYVMEKAGGMATTGKEAVLDVIPTDIHQRAPVILGSPDDVLEFLKVYEKHSQAQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	36.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001121100</u>
Locus ID:	2203


[View online »](#)

UniProt ID:	<u>P09467, Q2TU34</u>
RefSeq Size:	1546
Cytogenetics:	9q22.32
RefSeq ORF:	1014
Synonyms:	FBP
Summary:	Fructose-1,6-bisphosphatase 1, a gluconeogenesis regulatory enzyme, catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate. Fructose-1,6-diphosphatase deficiency is associated with hypoglycemia and metabolic acidosis. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Fructose and mannose metabolism, Glycolysis / Gluconeogenesis, Insulin signaling pathway, Metabolic pathways, Pentose phosphate pathway

Product images:



Coomassie blue staining of purified FBP1 protein (Cat# TP325497). The protein was produced from HEK293T cells transfected with FBP1 cDNA clone (Cat# [RC225497]) using MegaTran 2.0 (Cat# [TT210002]).