

Product datasheet for TP325547

Aldolase (ALDOA) (NM_001127617) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aldolase A, fructose-bisphosphate (ALDOA), transcript variant 4, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225547 protein sequence Red =Cloning site Green =Tags(s)

MPYQYPALTPEQKKELSDIAHRIVAPGKGILAADESTGSIKRLQSIGTENTEENRRFYRQLLLTADDRV
NPCIGGVILFHETLYQKADDGRFPQVIKSKGGVVGKVDKGWPLAGTNGETTTQGLDGLSERCAQYKK
DGADFAKWRCVLKIGEHTPSALAIMENANVLARYASICQONGIVPIVEPEILPDGDHDLKRCQYVTEKVL
AAVYKALSDHHIYLEGTLLKPNMVTPGHACTQKFSHEEIAMATVTALRRTPPPAVTGITFLSGGQSEEEA
SINLNAINKCPLLKPWALTFSYGRALQASALKAWGGKKENLAAQEEYVKRALANSLACQGKYTPSGQAG
AAASESLFVSNHAY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	39.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001121089</u>



[View online »](#)


Locus ID: 226
UniProt ID: [P04075](#)
RefSeq Size: 1630
Cytogenetics: 16p11.2
RefSeq ORF: 1092
Synonyms: ALDA; GSD12; HEL-S-87p

Summary: This gene encodes a member of the class I fructose-bisphosphate aldolase protein family. The encoded protein is a glycolytic enzyme that catalyzes the reversible conversion of fructose-1,6-bisphosphate to glyceraldehyde 3-phosphate and dihydroxyacetone phosphate. Three aldolase isozymes (A, B, and C), encoded by three different genes, are differentially expressed during development. Mutations in this gene have been associated with Glycogen Storage Disease XII, an autosomal recessive disorder associated with hemolytic anemia. Disruption of this gene also plays a role in the progression of multiple types of cancers. Related pseudogenes have been identified on chromosomes 3 and 10. [provided by RefSeq, Sep 2017]

Protein Families: Druggable Genome

Protein Pathways: Fructose and mannose metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Pentose phosphate pathway

Product images:



Coomassie blue staining of purified ALDOA protein (Cat# TP325547). The protein was produced from HEK293T cells transfected with ALDOA cDNA clone (Cat# [RC225547]) using MegaTran 2.0 (Cat# [TT210002]).