

Product datasheet for TP325705L

C13orf31 (LACC1) (NM_001128303) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 13 open reading frame 31 (C13orf31), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225705 protein sequence Red =Cloning site Green =Tags(s)

MAEAVLIDLFGLKLNQKNCHQTLKTLNAVQYHHAAKAKFLCIMCCSNISYERDGEQDNCEIETSNGLS
ALLEEFEIVSCPSMAATLYTIKQKIDEKNLSSIKVIVPRHRKTLMKAFIDQLFTDVYNFEFEDLQVTRG
GLFKQSIEINVITAQELRGIQNEIETFLRSLPALRGKLTITSSLIPDIFIHGFTTRTGGISYIPTLSSF
NLFSSSKRRDPKVVVQENLRRLANAAGFNVEKFYRIKTHHSNDIWIMGRKEPDSYDGITTNQRGVIAAL
GADCIPVIFADPVKKACGVAHAGWKGTLGVMATVNMIAEYGCSELDIVVWLGPSVGPCCFTLPRESA
EAFHNLHPACVQLFDSNPNCIDIRKATRILLEQGGILPQNIQDQNDLNLCTSCHPDKFFSHVRDGLNFG
TQIGFISIKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	47.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001121775](#)

Locus ID: 144811

UniProt ID: [Q8IV20](#)

RefSeq Size: 4288


Cytogenetics: 13q14.11

RefSeq ORF: 1290

Synonyms: C13orf31; FAMIN; JUVAR

Summary: This gene encodes an oxidoreductase that promotes fatty-acid oxidation, with concomitant inflammasome activation, mitochondrial and NADPH-oxidase-dependent reactive oxygen species production, and bactericidal activity of macrophages. The encoded protein forms a complex with fatty acid synthase on peroxisomes and is thought to be modulated by peroxisome proliferator-activated receptor signaling events. Naturally occurring mutations in this gene are associated with inflammatory bowel disease, Behcet's disease, leprosy, ulcerative colitis, early-onset Crohn's disease, and systemic juvenile idiopathic arthritis. [provided by RefSeq, Apr 2017]

Product images:



Coomassie blue staining of purified LACC1 protein (Cat# [TP325705]). The protein was produced from HEK293T cells transfected with LACC1 cDNA clone (Cat# [RC225705]) using MegaTran 2.0 (Cat# [TT210002]).