

Product datasheet for **TP325784L**

LSM14A (NM_001114093) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human LSM14A, SCD6 homolog A (<i>S. cerevisiae</i>) (LSM14A), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225784 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSGGTPYIGSKISLISKAEIRYEGILYTIDTENSTVALAKVRSFGTEDRPTDRPIPPRDEVFEYIIFRGS DIKDLTVCEPPKPQCSLPQDPAIVQSSLGSSTSSFQSMGSYGPFGRMPTYSQFSPSSLVGGQFGAVGVAG SSLTSFGTETSNSGTLQSSAVGSAFTQDTRSLKTQLSQGRSSPQLDPLRKSPTEQAVQTASAHLPA AVGRRSPVSTRPLPSASQKAGENQEHRRAEVHKVSRPENEQLRNDNKRQVAPGAPSAPRRGRGGHRRG RG RFGIRRDGPMKFEKDFDFESANAQFNKEEIDREFHNKLKLEDKLEKQEKPVNGEDKGD SGVDTQNSEG N ADEEDPLGPNCYYDKTKSFFDNISCDNRERRPTWAEERRLNAETFGIPLRPNRGRGGYRGRGGLGFRGG RGRGGGRGGTFTAPRGFRGGFRGGRGGRGGREFADFEYRKDNKVAA</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	50.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001107565</u>
Locus ID:	26065
UniProt ID:	<u>Q8ND56</u>
RefSeq Size:	3772
Cytogenetics:	19q13.11
RefSeq ORF:	1389
Synonyms:	C19orf13; FAM61A; RAP55; RAP55A
Summary:	Sm-like proteins were identified in a variety of organisms based on sequence homology with the Sm protein family (see SNRPD2; 601061). Sm-like proteins contain the Sm sequence motif, which consists of 2 regions separated by a linker of variable length that folds as a loop. The Sm-like proteins are thought to form a stable heteromer present in tri-snRNP particles, which are important for pre-mRNA splicing.[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified LSM14A protein (Cat# [TP325784]). The protein was produced from HEK293T cells transfected with LSM14A cDNA clone (Cat# [RC225784]) using MegaTran 2.0 (Cat# [TT210002]).