

Product datasheet for TP325796

ABLIM2 (NM_001130088) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human actin binding LIM protein family, member 2 (ABLIM2), transcript variant 7, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225796 representing NM_001130088 Red =Cloning site Green =Tags(s)
	MSAVSQQAAPSPLEKSPSTAILCNTCGNVCKGEVLRVQDKYFHIKCFVCKACGCDLAEGGFFVRQGEYI CTLDYQRLYGTRCFSCDQFIEGEVVSALGKTYHPDCFVCAVCRLPFPPGDRVTFNGKECMCQKCSLPVSV GSSAHLSQLRSCGGCGTEIKNGQALVALDKHWHLGCFKCKSCGKLLNAEYISKDGLPYCEADYHAKFGI RCDSCKEYITGRVLEAGEKHYPSCALCVRCQMFAEGEEMYLQGSSIWHPACRQAARTEDRNKETRTSS ESIISVPASSTSGSPSRVIYAKLGGEILDYRDLAALPKSKAIYDIDRPDMISYSPYISHSAGDRQSYGES PQLLSPTPTEGDQDDRSYQCRTSSPSSTGSVSLGRYTPTSRSPQHYSRPAARRSDGEDGSLDQDNRKQK SSWLMLKGDADTRTNSPDLDTQSLSHSSGTDRDPLQRMAGDSFHSREWFF
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	51.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001123560](#)

Locus ID: 84448

UniProt ID: [Q6H8Q1](#)


Cytogenetics: 4p16.1

RefSeq ORF: 1410

Summary: May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.[UniProtKB/Swiss-Prot Function]

Protein Pathways: Axon guidance

Product images:



Coomassie blue staining of purified ABLIM2 protein (Cat# TP325796). The protein was produced from HEK293T cells transfected with ABLIM2 cDNA clone (Cat# [RC225796]) using MegaTran 2.0 (Cat# [TT210002]).