

## Product datasheet for **TP325800**

### CD36 (NM\_001127444) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens CD36 molecule (thrombospondin receptor) (CD36), transcript variant 5, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225800 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MGCDRNCGLIAGAVIGAVLAVFGGILMPVGDLLIQKTIKKQVWLEEGTIAFKNWWKTGTEVYRQFWIFDV  
QNPQEVMMNSSNIQVKQRGPYTYRVRFLAKENVTQDAEDNTVSFLQPNGAIFEPSLSVGTEADNFTVLNL  
AVAAASHIYQNQFVQMILNSLINKSSMFQVRTLRELLWGYRDPFLSLVPYPVTTTVGLFPYNNTADG  
VYKVFNGKDNISKVAIIDTYKGRNLSYWESHCDMINGTDAASFPPFVEKSQVLQFFSSDICRSIYAVFE  
SDVNLKGIPVYRFLPSKAFASPVENPDNYCFCTEKIISKNCTSYGVLDISKCKEGRPVYISLPHFLYAS  
PDVSEPIDGLNPNEEEHRTYLDIEPITGFTLQFAKRLQVNLKPKSEKIQLKLNKRNIVIPILWLNETG  
TIGDEKANMFRSQVTGKINLLGLIEMILLSVGVMMFVAFMISYCACRSKTIK

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	52.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_001120916](#)

Locus ID: 948

UniProt ID: [P16671](#), [A4D1B1](#)

RefSeq Size: 1989

Cytogenetics: 7q21.11

RefSeq ORF: 1416


Synonyms: BDPLT10; CHDS7; FAT; GP3B; GP4; GPIV; PASIV; SCARB3

**Summary:** The protein encoded by this gene is the fourth major glycoprotein of the platelet surface and serves as a receptor for thrombospondin in platelets and various cell lines. Since thrombospondins are widely distributed proteins involved in a variety of adhesive processes, this protein may have important functions as a cell adhesion molecule. It binds to collagen, thrombospondin, anionic phospholipids and oxidized LDL. It directly mediates cytoadherence of Plasmodium falciparum parasitized erythrocytes and it binds long chain fatty acids and may function in the transport and/or as a regulator of fatty acid transport. Mutations in this gene cause platelet glycoprotein deficiency. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2014]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Adipocytokine signaling pathway, ECM-receptor interaction, Hematopoietic cell lineage, PPAR signaling pathway

### Product images:



Coomassie blue staining of purified CD36 protein (Cat# TP325800). The protein was produced from HEK293T cells transfected with CD36 cDNA clone (Cat# [RC225800]) using MegaTran 2.0 (Cat# [TT210002]).