

Product datasheet for **TP325878M**

Nuclear Factor 1 (NFIA) (NM_001134673) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nuclear factor I/A (NFIA), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225878 representing NM_001134673 Red=Cloning site Green=Tags(s)

MYSPLCLTQDEFHPFIEALLPHVRAFAYTWFNLQARKRKYFKKHEKRMSKEEERAVKDELLSEKPEVKQK
WASRLAKLRKDIRPEYREDFVLTVTGKKPPCCVLSNPDQKGMRRIDCLRQADKVVRLDLVMVILFKGI
PLESTDGERLVKSPQCSNPGLCVQPHHIGVSVKELDLYLAYFVHAADSSQSESPSQPSDADIKDQPENGH
LGFQDSFVTSVFSVTELVRSQTPIAAGTGPNFSLSDLESSSYSMSPGAMRRSLPSTSSTSKRLKS
VEDEMDSPEEPFYTGQGRSPGSGSQSSGWHEVEPGMPSPTTLKKSEKSGFSSPSPSQTSSLGTAFTQHH
RPVITGPRASPHATPSTLHFPTSPHQQPGPYFHPAIRYHPQETLKEFVQLVCPDAGQQAGQVGFNLNPN
GSSQGGKVNHPFLPTPMLPPPPPPPMARVPLPVPDTPKPTTSTEGGAASPTSPYSTPSTSPANRFVSVG
PRDPSFVNIPQQTQSWYLG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	55.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001128145](#)

Locus ID: 4774

UniProt ID: [Q12857](#)

Cytogenetics: 1p31.3


RefSeq ORF: 1527

Synonyms: BRMUTD; CTF; NF-I/A; NF1-A; NFI-A; NFI-L

Summary: This gene encodes a member of the NF1 (nuclear factor 1) family of transcription factors. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified NFIA protein (Cat# [TP325878]). The protein was produced from HEK293T cells transfected with NFIA cDNA clone (Cat# [RC225878]) using MegaTran 2.0 (Cat# [TT210002]).