

Product datasheet for TP325905L

CCDC149 (NM_001130726) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 149 (CCDC149), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225905 representing NM_001130726 Red =Cloning site Green =Tags(s)

MEEEAMNGDRTESDWQGLVSEYLVCKRKLESKKEALLILSKELDTCCQERDQYKLMANQLRERHQSLKK
 K
 YRELIDGDPSPPEKQANLAQLLRDSQDRNKHLEEKELQQLRGEVQGDNKLRLMTIAKQRLGDEAI
 GVRHFAAHEREDLVQQLERAKEQIESLEHDLQASVDELQDVKEERSSYQDKVERLNQELNHILSGHENRI
 IDVDALCMENRYLQERLQQLHEEVNLLKSNIKYKNALERRKNSKGQKSSSSALTGVLSAKQVQDLLSE
 DHGCSLPATPQSISDLKSLATALLEIHEKNMVIQHQRQTNKILGNRVAELEKKLRTLEVSGLWVSLPGGK
 DTILFSDPTLPSGQRSRSPLLKFVEQPTENKADPKDGEAQKQEEDESCAAAEALTAPEDAGRPVANSAN
 QSRGNQCKLFHPSLPQLPSEEEVNSLGREIKLTKEQAAAELEEVRRSPIEGQRSETGPAPPGLAIQGE
 LPKSHLDSFEASRPAKASTPEDGKGIPGEGGMRSTVKT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	59.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001124198</u>
Locus ID:	91050
UniProt ID:	<u>Q6ZUS6</u>
Cytogenetics:	4p15.2
RefSeq ORF:	1587

Product images:



Coomassie blue staining of purified CCDC149 protein (Cat# [TP325905]). The protein was produced from HEK293T cells transfected with CCDC149 cDNA clone (Cat# [RC225905]) using MegaTran 2.0 (Cat# [TT210002]).