

Product datasheet for **TP325961**

NT5C2 (NM_001134373) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human 5'-nucleotidase, cytosolic II (NT5C2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC225961 protein sequence Red=Cloning site Green=Tags(s)

MSTSWSDRLQNAADMPANMDKHALKKYRREAYHRVFNRLAMEKIKCFGDMDYTLAVYKSPEYESLGF
ELTVERLVSIGYPQELLSFAYDSTFPTRGLVFDTLYGNNLLKVDAYGNLLVCAHGFNFIRGPETREQYPNK
FIQRDDTERFYILNTLFLNPETYLLACLVDFFTNCPRYTSCETGFKDGDLMFSYRSMFQDVRDAVDWVHY
KGSLKEKTVENLEKYVVKDGKLPDLLSRMKEVGKVFATNSDYKYTDKIMTYLFDFFHGPKPGSSHRPWQ
SYFDLILVDARKPLFFGEGTVLRQVDTKTGKLIKITYTGPLQHGVYSGGSSDTICDLLGAKGKDILYIG
DHIFGDILKSKKRQGWRTFLVIPELAQELHVWTDKSSLFEELQSLDIFLAELYKHLDSSSNERPDISSIQ
RRIKKVTHDMDMCMGMMGSLFRSGSRQTLFASQVMRYADLYAASFINLLYPFSYLFRAAHVLMPHSTV
EHTHVDINEMESPLATRNRTSVDFKDTDYKRHLQRSISEIKPPNLFPLAPQEITHCHDEDDDEEEEEEE
E

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	64.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001127845
Locus ID:	22978
UniProt ID:	P49902 , A0A384MED8 , A8K6K2
RefSeq Size:	3453
Cytogenetics:	10q24.32-q24.33
RefSeq ORF:	1686
Synonyms:	cN-II; GMP; NT5B; PNT5; SPG45; SPG65
Summary:	This gene encodes a hydrolase that serves as an important role in cellular purine metabolism by acting primarily on inosine 5'-monophosphate and other purine nucleotides. [provided by RefSeq, Oct 2011]
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism, Purine metabolism, Pyrimidine metabolism

Product images:



Coomassie blue staining of purified NT5C2 protein (Cat# TP325961). The protein was produced from HEK293T cells transfected with NT5C2 cDNA clone (Cat# [RC225961]) using MegaTran 2.0 (Cat# [TT210002]).