

## Product datasheet for TP326168

### Breast cancer suppressor candidate 1 (VWA5A) (NM\_001130142) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human von Willebrand factor A domain containing 5A (VWA5A), transcript variant 3, 20 µg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >RC226168 representing NM\_001130142  
**Red**=Cloning site **Green**=Tags(s)

MVHFCGLLTLHREPVLKSSISVSVNIYEFVAGVSATLNYENEEKVPLEAFFVFPMDEDSAVYSFEALVDG  
KKIVAEQLQDKMKARTNYEKAISQGHQAFLLLEGDSSSRDVFSCNVGNLQPGSKAAVTLKYVQELPLEADGA  
LRFVLPVAVLNPRYQFSGSSKDSCLNVKTPIVPVEDLPYTLMSVATIDSQHGIEKVQSNCLPSPEYLGED  
KTSAQVSLAAGHKFDRDVELLIYYNEVHTPSVVLEMGPMPNPKPGHLMGDPMSAMVSFYPNIPEDQPSNTCG  
EFIFLMDRSGSMQSPMSSQDTSQRLRIQAAKETLILLKSLPIGCFNYIYFGSSYEACFPESVKYTQQT  
EEALGRVKLMQADLGGTEILAPLQNIYRGPSIPGHPLQLFVFTDGEVTDTFSVIKEVRINRQKHRCFSFG  
IGEGTSTSLIKGIARASGGTSEFITGKDRMQSKALRTLKRSLQPVVEDVSLSWHLPPGLSAKMLSPEQTV  
IFRGQRLISYAQLTGRMPAAETTGEVCLKYTLQGKTFEDKVTFLPQPKPDVNLTIHRLAAKSLQTKDMG  
LRETPASDKKDALNLSLESGVISSFTAFIANKELNKPVQGPLAHRDVPRPILLGASAPLKIKCQSGFRK  
ALHSDRPPSASQPRGELMCYKAKTFQMDDYSLCGLISHKDQHSPGFGENHLVQLIYHQNANGSWDLNEDL  
AKILGMSLEEIMAAQPAELVDSSGWATILAVIWLHSNGKDLKCEWELLERKAVAWMRAHAGSTMPSSVKA  
AITFLKSSVDP AIFAF

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK

**Predicted MW:** 86.3 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining


**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001123614</a>
<b>Locus ID:</b>	4013
<b>UniProt ID:</b>	<a href="#">Q00534</a> , <a href="#">A0A024R3H3</a>
<b>Cytogenetics:</b>	11q24.2
<b>RefSeq ORF:</b>	2358
<b>Synonyms:</b>	BCSC-1; BCSC1; LOH11CR2A
<b>Summary:</b>	May play a role in tumorigenesis as a tumor suppressor. Altered expression of this protein and disruption of the molecular pathway it is involved in, may contribute directly to or modify tumorigenesis.[UniProtKB/Swiss-Prot Function]

**Product images:**

Coomassie blue staining of purified VWA5A protein (Cat# TP326168). The protein was produced from HEK293T cells transfected with VWA5A cDNA clone (Cat# [RC226168]) using MegaTran 2.0 (Cat# [TT210002]).