

Product datasheet for **TP326414M**

DNMT1 (NM_001130823) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DNA (cytosine-5-)-methyltransferase 1 (DNMT1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T



[View online »](#)

Expression cDNA >RC226414 representing NM_001130823
Clone or AA **Red**=Cloning site **Green**=Tags(s)
Sequence:



MPARTAPARVPTLAVPAISLPDDVRRRLKDLERDSLTEKECVKEKLNLLHEFLQTEIKNQLCDLETKLK
 EELSEEGYLAKVKSLLNKDLLENGAHAYNREVNGLRENGNQARSEARRVGMADANSPPKPLSKPRTPRR
 SKSDGEAKRSRDPASASQVTGIRAEPSPRITRKSTRQTTITSHFAKGPAAKRPQEESERAKSDESIK
 EEDKDQDEKRRRVTSRERVARPLPAEPPERAKSGTRTEKEEERDEKEEKRLRSQTKPTPKQKLKEEPDR
 EARAGVQADEDGDEKDEKKHRSQPKDLAAKRRPEEKEPEKVNQPISDEKDEDEKKEEKRRKTPKEPTE
 KKMARAKTVMNSKTHPPKCIQCGQYLDDPDLKYQHPPDAVDEPQMLTNEKLSIFDANESGFESYEALPQ
 HKLTCFSVYCKHGHLCPIDTGLIEKNIELFFSGSAKIYDDPSLEGGVNGKNLGPINEWWTGFDGGGK
 ALIGFSTFAEYILMDPSPEYAPIFGLMQEKIYISKIVVEFLQSNSDSTYEDLINKIETTVPSSGLNLNR
 FTEDSLLRHAQFVVEQVESYDEAGDSDEQPIFLTPCMRDLIKLAGVTLGQRRRAQARRQTIHSTREKDRG
 PTKATTTKLQYQIFDTFFAEQIEKDDREDKENAFKRRRCGVCVCQQPECGKCKACKDMVKFGGSGRSKQ
 ACQERRCPNMAMKEADDDDEVDNIPEMPSPKMHQGKKKKQKNKRISWVGEAVKTDGKKSYYKKVCIDA
 ETLEVGDVSVIPDDSSKPLYLARVTALWEDSSNGQMFHAWFCAGTDTVLGATSDPLELFLVDECEDMQ
 LSYIHSKVKVIYKAPSENWAMEGGMDPESLLEGDDGKTYFYQLWYDQDYARFESPPKTQPTEDNKFKFCV
 SCARLAEMRQKEIPRVLEQLEDLSRVLYSATKNGILYRVGDGVYLPPEAFTFNILSSPVKRPKEPV
 DEDLYPEHYRKYSDYIKGSNLDAPEPYRIGRIKEIFCPKKSNGRPNETDIKIRVNKFYRPENTHKSTPAS
 YHADINLLYWSDEEAVVDFKAVQGRCTVEYGEDLPECVQVYSMGGPNRFYFLEAYNAKSKSFEDPPNHAR
 SPGNKKGKGGKGGKPKSQACEPSEPEIEIKLPKLRDLDFVSGCGGLSEGFHQAGISDTLWAIEMWDPAA
 QAFRLNPNPGSTVFTEDCNILLKLVMAGETTNSRGQRLPQKGDVEMLCGGPPCQGFSGMNRFNRSRTYSKFK
 NSLVVSFLSYCDYYRPRFFLEENVRNFVSFKRSMVLKLTLRCLVRMGYQCTFGVLQAGQYGAQTRRRAI
 ILAAAPGEKLPFPEPLHVFAPRACQLSVVDDKKFVSNITRLSSGPFRTITVRDMSDLPEVRNGASAL
 EISYNGEPQSWFQRQLRGAQYQPILRDHICKDMSALVAARMRHIPLAPGSDWRDLNIEVRLSDGTMARK
 LRYTHHDRKNGRSSSGALRGVCSCVEAGKACDPAARQFNTLIPWCLPHTGNRHNHWAGLYGRLEWDGFFS
 TTVTNPEPMGKQGRVLHPEQHRVSVRECARSQGFPDTRFLFGNILDKHRQVGNVAVPPPLAKAIGLEIKL
 CMLAKARESASAKIKEEEAAKD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 184.6 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Bioactivity:	DNMT1 activity verified in a biochemical assay: DNMT1 (DNA (cytosine-5-)-methyltransferase 1) (TP326414) is a key methyltransferase that is responsible for maintaining methylation patterns established in development. DNMT1 preferentially methylates hemi-methylated CpG dinucleotides and associates with DNA replication sites in S phase maintaining the methylation pattern in the newly synthesized strand. Varying concentrations of DNMT1 were added to a microplate containing a bound methyltransferase substrate. After incubation, the resulting methylated DNA residues were detected immunologically and a colorimetric signal was generated and measured.
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001124295
Locus ID:	1786
UniProt ID:	P26358 , I6L9H2
Cytogenetics:	19p13.2
RefSeq ORF:	4896
Synonyms:	ADCADN; AIM; CXXC9; DNMT; HSN1E; m.HsaI; MCMT
Summary:	This gene encodes an enzyme that transfers methyl groups to cytosine nucleotides of genomic DNA. This protein is the major enzyme responsible for maintaining methylation patterns following DNA replication and shows a preference for hemi-methylated DNA. Methylation of DNA is an important component of mammalian epigenetic gene regulation. Aberrant methylation patterns are found in human tumors and associated with developmental abnormalities. Variation in this gene has been associated with cerebellar ataxia, deafness, and narcolepsy, and neuropathy, hereditary sensory, type IE. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified DNMT1 protein (Cat# [TP326414]). The protein was produced from HEK293T cells transfected with DNMT1 cDNA clone (Cat# [RC226414]) using MegaTran 2.0 (Cat# [TT210002]).