

Product datasheet for **TP326508**

RAB34 (NM_001144943) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens RAB34, member RAS oncogene family (RAB34), transcript variant 8, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226508 representing NM_001144943 Red =Cloning site Green =Tags(s)

MSHLPGLELRREAPLLGPLLSPFPLPAGSWHRQMLRSSLRFITNSAGAPCKAAGRMINILAPVRRDRVL
AELPQCLRKEAALHGHKDFHPRVTCACQEHRTGTVGFKISKVIVVGDLSVGKTCLINRFCKDTDFDKNYKA
TIGVDFEMERFEVLGIPFSLQLWDTAGQERFKCIASTYYRGAQAIIVFNLNDVASLEHTKQWLADALKE
NDPSSVLLFLVGSKKDLSTPAQYALMEKDALQVAQEMKAQYWAQVSSLTGENVREFFFRVAALTFEANVLA
ELEKSGARRIGDVVRINSDDSNLYLTASKKKPTCCP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	35 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138415
Locus ID:	83871



[View online »](#)

UniProt ID: [B4DNC0](#)

Cytogenetics: 17q11.2


RefSeq ORF: 948

Synonyms: NARR; RAB39; RAH

Summary: This gene encodes a protein belonging to the RAB family of proteins, which are small GTPases involved in protein transport. This family member is a Golgi-bound member of the secretory pathway that is involved in the repositioning of lysosomes and the activation of macropinocytosis. Alternative splicing of this gene results in multiple transcript variants. An alternatively spliced transcript variant produces the nine-amino acid residue-repeats (NARR) protein, which is a functionally distinct nucleolar protein resulting from a different reading frame. [provided by RefSeq, Dec 2016]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified RAB34 protein (Cat# TP326508). The protein was produced from HEK293T cells transfected with RAB34 cDNA clone (Cat# [RC226508]) using MegaTran 2.0 (Cat# [TT210002]).