

## Product datasheet for **TP326633M**

### **BAP31 (BCAP31) (NM\_001139441) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human B-cell receptor-associated protein 31 (BCAP31), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226633 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MSLQWTAVATFLYAEVFWLLLCIPFISPKRWQKIFKSRLVELLVSYGNTFFWLVILVLLVIDAVREI REYDDVTEKYNLQNNPGAMEHFHMKLFRAQRNLYIAGFSLLSFLLRRLVTLISQQATLLASNEAFKKQA ESASEAAKKYMEENDQLKKGAAVDGGKLDVGNAEVKLEENRSLKADLQKLKDELASTKQKLEKAENQVL AMRKQSEGLTKEYDRLLLEEHAQLAAVDGPMDDKKEE</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	27.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001132913</a></u>
Locus ID:	10134




[View online »](#)

UniProt ID: [P51572](#)  
RefSeq Size: 1346  
Cytogenetics: Xq28  
RefSeq ORF: 738  
Synonyms: 6C6-AG; BAP31; CDM; DDCH; DXS1357E

**Summary:** This gene encodes a member of the B-cell receptor associated protein 31 superfamily. The encoded protein is a multi-pass transmembrane protein of the endoplasmic reticulum that is involved in the anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi and in caspase 8-mediated apoptosis. Microdeletions in this gene are associated with contiguous ABCD1/DXS1375E deletion syndrome (CADD5), a neonatal disorder. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 16. [provided by RefSeq, Jan 2012]

**Protein Families:** Druggable Genome, Transmembrane

### Product images:



Coomassie blue staining of purified BCAP31 protein (Cat# [TP326633]). The protein was produced from HEK293T cells transfected with BCAP31 cDNA clone (Cat# [RC226633]) using MegaTran 2.0 (Cat# [TT210002]).