

## Product datasheet for **TP326640M**

### **MBD3L5 (NM\_001136507) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human methyl-CpG-binding domain protein 3-like 5-like (MBD3L5), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226640 representing NM_001136507 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MGEPAFTSFPSPVVLGKLRNMMPWALQKKREIHMakahrrraarsalpmrltscifrrpvtrirshpdn QVRRRKGDEHLEKPPQLCAYRRLQALQPCSSQGEgssplhlesvlsilapgtagesldragAERVRIPL PTPGRFPAVAGGPTPGMGCQLPPPLSGQLVTPADIRRQARRVKKARERLAKALQADRLARQAEMLTGG</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	22.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001129979</a></u>
Locus ID:	284428
UniProt ID:	<u><a href="#">A6NJ08</a></u>




[View online »](#)

Cytogenetics: 19p13.2

RefSeq ORF: 624

### Product images:



Coomassie blue staining of purified MBD3L5 protein (Cat# [TP326640]). The protein was produced from HEK293T cells transfected with MBD3L5 cDNA clone (Cat# [RC226640]) using MegaTran 2.0 (Cat# [TT210002]).