

Product datasheet for TP326656M

DSN1 (NM_001145318) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens DSN1, MIND kinetochore complex component, homolog (<i>S. cerevisiae</i>) (DSN1), transcript variant 5, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226656 representing NM_001145318 Red=Cloning site Green=Tags(s)

MSKTHDQHQLLESSLPVEVFAKTSASLEMNQGVSEERIHGSSPKKGGNCDLSHQERLQSKSLHLSPQEQS
ASYQDRRQSWRRASMKETNRRKSLHPIHQGITELSRISVDLAESKRLGCLLLSSFQFSIQKLEPFLRDT
KGFLESFRKASSLSEELKHFADGLETDGTLQKCFEDSNGKASDFSLEASVAEMKEYITKFSLERQTWD
QLLLHYQQEAKEILSRGSTEAKITEVKVEPMTYLGSSQNEVLNTPDYQKILQNQSKVFDCELMVMDLQ
GSVKQLQAFMDESTQCFQKVSQVQLGKRSMQQLDPSARKLLKLQLNPPAIHGSGSGSCQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	38.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138790
Locus ID:	79980



[View online »](#)

UniProt ID:	Q9H410
Cytogenetics:	20q11.23
RefSeq ORF:	1020
Synonyms:	C20orf172; dj469A13.2; hKNL-3; KNL3; MIS13
Summary:	This gene encodes a kinetochore protein that functions as part of the minichromosome instability-12 centromere complex. The encoded protein is required for proper kinetochore assembly and progression through the cell cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009]

Product images:



Coomassie blue staining of purified DSN1 protein (Cat# [TP326656]). The protein was produced from HEK293T cells transfected with DSN1 cDNA clone (Cat# [RC226656]) using MegaTran 2.0 (Cat# [TT210002]).