

## Product datasheet for TP326668

### HOMER3 (NM\_001145722) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human homer homolog 3 (Drosophila) (HOMER3), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226668 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MSTAREQPIFSTRAHVFQIDPATKRNWIPAGKHALTVSYFYDATRNVYRIISIGGAKAIINSTVTPNMTF  
TKTSQKFGQWADSRANTVYGLGFASEQHLTQFAEKFQEVKEAARLAREKSQDGGELTSPALGLASHQVPP  
SPLVSANGPGEEKLFRSQSADAPGPTERERLKKMLSEGSVGEVQWEAEFFALQDSNNKLAGALREANAAA  
AQWRQQLAQRAEAERLRQRVAELEQAASEVTPTTGEKEGLGQGSLEQLEALVQTKDQEIQTLKSQTGG  
PREALEAAEREETQQKVQDLETRNAELEHQLRAMERSLEEARAERERARAIEVGRAAQLLDVRLFELSELR  
EGLARLAEAAP

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	39.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001139194</a>




[View online »](#)

Locus ID: 9454  
UniProt ID: [Q9NSC5](#)  
RefSeq Size: 1859  
Cytogenetics: 19p13.11  
RefSeq ORF: 1083  
Synonyms: HOMER-3; VESL3

**Summary:** This gene encodes a member of the HOMER family of postsynaptic density scaffolding proteins that share a similar domain structure consisting of an N-terminal Enabled/vasodilator-stimulated phosphoprotein homology 1 domain which mediates protein-protein interactions, and a carboxy-terminal coiled-coil domain and two leucine zipper motifs that are involved in self-oligomerization. The encoded protein binds numerous other proteins including group I metabotropic glutamate receptors, inositol 1,4,5-trisphosphate receptors and amyloid precursor proteins and has been implicated in diverse biological functions such as neuronal signaling, T-cell activation and trafficking of amyloid beta peptides. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Mar 2009]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified HOMER3 protein (Cat# TP326668). The protein was produced from HEK293T cells transfected with HOMER3 cDNA clone (Cat# [RC226668]) using MegaTran 2.0 (Cat# [TT210002]).