

Product datasheet for **TP326700M**

EVPLL (NM_001145127) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human envoplakin-like (EVPLL), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC226700 representing NM_001145127
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MQASADQVERDILETQKRLQQDRLNSEQSQALQHQQETGSSLKEAEVLLKDLFLDVKARRLKHPQAEET
EKDIEQLHERVTQCEAIEYCALYKMLVPPRRGIQGRGTRAGAETEAGLRRPVWAGHGGAGGTDGRGAQHR
AEGDQRPRRAAAEPGGAGCRHHPEPIPRPTEGGVARAEPGQPVHALQGCTWQLSALAEQQRRLQQDWS
DLMADPAGVRREYEHFKQHELLSQEQSVNQLLEEDGKRMVELRHPAVGPIQAHQEALKMEWQNFLNLCICQ
ETQLQHVEDYSRILCPSSSPH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	33.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138599
Locus ID:	645027



[View online »](#)

UniProt ID: [A8MZ36](#)
Cytogenetics: 17p11.2
RefSeq ORF: 903

Product images:



Coomassie blue staining of purified EVPLL protein (Cat# [TP326700]). The protein was produced from HEK293T cells transfected with EVPLL cDNA clone (Cat# [RC226700]) using MegaTran 2.0 (Cat# [TT210002]).