

Product datasheet for TP326704L

GNDF Receptor alpha 1 (GFRA1) (NM_001145453) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GDNF family receptor alpha 1 (GFRA1), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226704 representing NM_001145453 Red=Cloning site Green=Tags(s)

MFLATLYFALPLDLLLSAEVSGGDRLCDVKASDQCLKEQSCSTKYRTLRLRQCVAGKETNFSLASGLEAKD
 ECRSAMEALKQKSLYNCRCKRGMKKEKNCLRIYWSMYQSLQGNDLLEDSPYEPVNSRLSDIFRVVPFISV
 EHIPKGNCLDAKACNLDDICKKYRSAYITPCTTSVSNVDCNRRKCHKALRQFFDKVPAKHSYGMLFCS
 CRDIACTERRRQTIVPVCSYEEREKPNCLNLQDSCKTNYICRSRLADFFTNCPESRSVSSCLKENYADC
 LLAYSGLIGTVMTPNYIDSSLSVAPWCDCSNSGNDLEECLKFLNFFKDNTCLKNAIQAFNGSDVTWVQ
 PAFPVQTTTATTTTALRVKNKPLGPAGSENIPTHVLPSCANLQAQKLKSNVSGNTHLCISNGNYEKEGL
 GASSHITTKSMAAPPSCGLSPLLVLVTALSTLLSLTETS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	48.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001138925](#)

Locus ID: 2674

UniProt ID: [P56159](#)

Cytogenetics: 10q25.3


RefSeq ORF: 1380

Synonyms: GDNFR; GDNFRA; GFR-ALPHA-1; GFRalpha-1; RET1L; RETL1; TRNR1

Summary: This gene encodes a member of the glial cell line-derived neurotrophic factor receptor (GDNFR) family of proteins. The encoded preproprotein is proteolytically processed to generate the mature receptor. Glial cell line-derived neurotrophic factor (GDNF) and neurturin (NTN) are two structurally related, potent neurotrophic factors that play key roles in the control of neuron survival and differentiation. This receptor is a glycosylphosphatidylinositol (GPI)-linked cell surface receptor for both GDNF and NTN, and mediates activation of the RET tyrosine kinase receptor. This gene is a candidate gene for Hirschsprung disease. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed. [provided by RefSeq, Jan 2016]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified GFRA1 protein (Cat# [TP326704]). The protein was produced from HEK293T cells transfected with GFRA1 cDNA clone (Cat# [RC226704]) using MegaTran 2.0 (Cat# [TT210002]).