

Product datasheet for TP326792L

C5orf24 (NM_001135586) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 5 open reading frame 24 (C5orf24), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226792 protein sequence Red =Cloning site Green =Tags(s) MMHPVASSNPAFCGPGKPSCLNEDAMRAADQFDIYSSQSKYSHTVNHKPMVCQRQDPLNETHLQTTSGR SIEIKDELKKKKNLNRSGKRGRPSGTTKSAGYRTSTGRPLGTTKAAGFKTSPGRPLGTTKAAGYKVSPGR PPGSIKALSRLADLGYGCGTAAFPYPMMHGRAVHGVEETSSEVKPPNE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129058
Locus ID:	134553
UniProt ID:	Q7Z6I8




[View online »](#)

RefSeq Size: 5083

Cytogenetics: 5q31.1

RefSeq ORF: 564

Product images:



Coomassie blue staining of purified C5orf24 protein (Cat# [TP326792]). The protein was produced from HEK293T cells transfected with C5orf24 cDNA clone (Cat# [RC226792]) using MegaTran 2.0 (Cat# [TT210002]).