

Product datasheet for **TP326813L**

ICA1 (NM_001136020) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens islet cell autoantigen 1, 69kDa (ICA1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226813 representing NM_001136020 Red =Cloning site Green =Tags(s)

MSGHKCYPWDLQDRYAQDKSVVNKMQQKYWETKQAFIKATGKKEDHWASDADLDAKLELFHSIQRTCL
DLSKAIVLYQKRICFLSQEENELGKFLRSQGFQDKTRAGKMMQATGKALCFSSQQLALRNPLCRFHQEV
ETFRHRAISDTWLTVNRMEQCRTEYRGALLWMKDVQSQELDPDLYKQMEKFRKVQTVRLAKKNFDKMKMD
VCQKVDLLGASRCNLLSHMLATYQTLLHFWEKTSHTMAAIHESFKGYQPYEFTTLKSLQDPMKKLVEKE
EKKKINQQESTDAAVQEPSQLISLEEENQRKESSFKTEDGKSILSALDKGSTHTACSGPIDELDMKSE
EGACLGPVAGTPEPEGADKDDLLLLSEIFNASSLEEGEFSKEWAAVFGDGQVKEPVPTMALGEPDPAQT
GSGFLPSQLLDQNMKDLQASLQEPAKAASDLTAWFSLFADLDPLSNPDAVKGKTDKEHELLNA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	54.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129492



[View online »](#)

Locus ID: 3382

UniProt ID: [Q05084](#), [A0A024RA29](#)

Cytogenetics: 7p21.3


RefSeq ORF: 1244

Synonyms: ICA69; ICAp69

Summary: This gene encodes a protein with an arfaptin homology domain that is found both in the cytosol and as membrane-bound form on the Golgi complex and immature secretory granules. This protein is believed to be an autoantigen in insulin-dependent diabetes mellitus and primary Sjogren's syndrome. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Feb 2013]

Protein Pathways: Type I diabetes mellitus

Product images:



Coomassie blue staining of purified ICA1 protein (Cat# [TP326813]). The protein was produced from HEK293T cells transfected with ICA1 cDNA clone (Cat# [RC226813]) using MegaTran 2.0 (Cat# [TT210002]).