

Product datasheet for TP326861M

GPR56 (ADGRG1) (NM_001145771) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human G protein-coupled receptor 56 (GPR56), transcript variant 4, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226861 protein sequence Red=Cloning site Green=Tags(s)

MTPQSLQTTFLLSLLFLVQGAHGRGHREDFRCSQRNQTHRSSLHYKPTPDLRISIENSEEALTVHAP
FPAAHPASRSFPDPRGLYHFCLYWNRHAGRLHLLYGKRDFLLSDKASSLLCFQHQEESLAQGPELLATS
TSWWSPQNIPLSAASFTFSFSPHTAAHNASVDMCELKRDLQLLSQFLKHPQKASRRPSAAPASQQLQ
SLESKLTSVRFMGDMVSFEEDRINATVWKLQPTAGLQDLHIHSRQEEEQSEIMEYSVLLPRTLQRTKGR
SGEAEKRLLLVDFSSQALFQDKNSSHVLGEKVLGIVVQNTKVANLTPWVLTQFHLQPKNVTLQCVFWV
EDPTLSSPGHWSSAGCETVRRETQTSCFCNHLTYFAVLMVSSVEVDVAVHKHYLSLLSYVGCVVSALACL
TIAAYLCSRVP LPCRKRPRDYTIKVMNLLLAVFLDTSFLLSEPVALTGSEAGCRASAI FLHFSLTCL
SWMGLEGNLYRLVVEVFGTYVPGYLLKLSAMGWGFPIFLVTLVALVDVDNYGPIILAVHRTPEGVIYPS
MCWIRDSLVSITNLGLFSLVFLNMAMLATMVVQILRLRPHTQKWSHVLTLLGLSLVLGLPWALIFFSF
ASGTFQLVVLVYLSIITSFQGLFIFIWYWSMRLQARGGPSPLKSNSDSARLPISSGSTSSRI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	75 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001139243
Locus ID:	9289
UniProt ID:	Q9Y653 , A0A024R6U7
RefSeq Size:	4269
Cytogenetics:	16q21
RefSeq ORF:	2079
Synonyms:	BFPP; BPPR; GPR56; TM7LN4; TM7XN1
Summary:	This gene encodes a member of the G protein-coupled receptor family and regulates brain cortical patterning. The encoded protein binds specifically to transglutaminase 2, a component of tissue and tumor stroma implicated as an inhibitor of tumor progression. Mutations in this gene are associated with a brain malformation known as bilateral frontoparietal polymicrogyria. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014]
Protein Families:	Druggable Genome, GPCR, Transmembrane

Product images:



Coomassie blue staining of purified ADGRG1 protein (Cat# [TP326861]). The protein was produced from HEK293T cells transfected with ADGRG1 cDNA clone (Cat# [RC226861]) using MegaTran 2.0 (Cat# [TT210002]).