

Product datasheet for **TP326864M**

LMO2 (NM_001142316) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human LIM domain only 2 (rhombotin-like 1) (LMO2), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226864 representing NM_001142316 Red =Cloning site Green =Tags(s) MSSAIERKSLDPSEEPVDEVLQIPPSLLTCGGCQQNIGDRYFLKAIDQYWHEDCLSCDLCGCRLEGEVGRRLYYKLGKRLCRRDYLRFLFGQDGLCASCDKRIRAYEMTMRVKDKVYHLECFKCAACQKHFCVGDYRLLINSDIVCEQDIYEWTKINGMI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001135788</u>
Locus ID:	4005
UniProt ID:	<u>P25791</u>



[View online »](#)

Cytogenetics:	11p13
RefSeq ORF:	474
Synonyms:	LMO-2; RBTN2; RBTNL1; RHOM2; TTG2
Summary:	LMO2 encodes a cysteine-rich, two LIM-domain protein that is required for yolk sac erythropoiesis. The LMO2 protein has a central and crucial role in hematopoietic development and is highly conserved. The LMO2 transcription start site is located approximately 25 kb downstream from the 11p13 T-cell translocation cluster (11p13 ttc), where a number T-cell acute lymphoblastic leukemia-specific translocations occur. Alternative splicing results in multiple transcript variants encoding different isoforms.[provided by RefSeq, Nov 2008]
Protein Families:	Druggable Genome

Product images:

Coomassie blue staining of purified LMO2 protein (Cat# [TP326864]). The protein was produced from HEK293T cells transfected with LMO2 cDNA clone (Cat# [RC226864]) using MegaTran 2.0 (Cat# [TT210002]).