

Product datasheet for **TP326878**

CXorf49B (NM_001145139) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome X open reading frame 49B (CXorf49B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226878 representing NM_001145139 Red =Cloning site Green =Tags(s) MSSPDKVSVCGAGFDLEGGKKAGSRTASPGAPGAHSHGLDLGVPGSGDGKSESGFTDPEGFSFESESELI EQGRVVLWGREGRPVDDQGDWVDYSFYLADEPAAIVPPPSVQGHPPFEGAAAEGSAENWADAIEVGP GRDVLGHSPGKWQQASAGRLHLCGPGVRAWKNPERGSKSRWSLRVDPQQPSAKGPTRLPHDSDSADES SDLPLMKVIGICRNEGSQAKPGSPKCRADTSRQASFHCKESYLPVPGRFLTSAPRGLTPVAERPAVGELED SPQKKMQSRAWGKVEVRPSCSGAAAAGALPQGLSRRKMAGGKSLGGASQLALGRGFACGERLSAAPPE PATFPPFSGVRPQGMSSKPKPKHSSPGKKPAGRKTRESQAAAREDNPNRDEVPRALPTHRPLRLS VRRGEFSSSDPNIRAPQLPGTSEPSAYSPGGLVPRRHAPSGNQPPVHPPRPERQQPPGAQGCPRCIWL QREIEDLTQQLAAMQFLTDKFDL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	54.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001138611](#)

Locus ID: 100132994


UniProt ID: [A8MYA2](#)

Cytogenetics: Xq13.1

RefSeq ORF: 1542

Synonyms: CXorf49

Product images:



Coomassie blue staining of purified CXorf49B protein (Cat# TP326878). The protein was produced from HEK293T cells transfected with CXorf49B cDNA clone (Cat# [RC226878]) using MegaTran 2.0 (Cat# [TT210002]).