

Product datasheet for TP326950L

CNGB1 (NM_001135639) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cyclic nucleotide gated channel beta 1 (CNGB1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC226950 representing NM_001135639 Red=Cloning site Green=Tags(s)
	<p>MLGWWQRVLPQPPGTPRKTKMQEEEEVEPEPEMEAEVEPEPNPEEAETES SMPPEESFKEEEVAVADPS PQETKEAALTSTISLRAQGAEISEMNSPSRRVLTWLMKGVKQVIPQPVHSITEDPAQILGHGSTGDTGCT DEPNEALEAQDTRPGLRLLLWLEQNLERVLPQPPKSSEVWRDEPAVATGAASDPAPPGRPQEMGPKLQAR ETPSLPTPIPLQPKEEPKEAPAPEPQPGSQAQTSSLPPTRDPARLVAVVLRHLEMALPQPVLHGKIGEQE PDSPGICDVQTRVMGAGGL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	32.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001129111</u>
Locus ID:	1258



[View online »](#)

UniProt ID: [Q14028](#)

Cytogenetics: 16q21

RefSeq ORF: 897


Synonyms: CNCG2; CNCG3L; CNCG4; CNG4; CNGB1B; GAR1; GARP; GARP2; RCNC2; RCNCb; RCNCbeta; RP45

Summary: In humans, the rod photoreceptor cGMP-gated cation channel helps regulate ion flow into the rod photoreceptor outer segment in response to light-induced alteration of the levels of intracellular cGMP. This channel consists of two subunits, alpha and beta, with the protein encoded by this gene representing the beta subunit. Defects in this gene are a cause of cause of retinitis pigmentosa type 45. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]

Protein Families: Druggable Genome, Ion Channels: Cyclic nucleotide gated

Protein Pathways: Olfactory transduction

Product images:



Coomassie blue staining of purified CNGB1 protein (Cat# [TP326950]). The protein was produced from HEK293T cells transfected with CNGB1 cDNA clone (Cat# [RC226950]) using MegaTran 2.0 (Cat# [TT210002]).