

## Product datasheet for **TP327004**


### **C4orf48 (NM\_001141936) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 4 open reading frame 48 (C4orf48), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227004 representing NM_001141936 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MAPPPACRSPMSPPPLLLLLLLSLALLGARARAEPAGSAVPAQSRPCVDCHAFEFMQRALQDLRKTACS LDARTETLLQAERRALCACPAGH  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	10 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001135408</a></u>
Locus ID:	401115
Cytogenetics:	4p16.3
RefSeq ORF:	285
Synonyms:	CHR4_55



[View online »](#)

**Product images:**


Coomassie blue staining of purified C4orf48 protein (Cat# TP327004). The protein was produced from HEK293T cells transfected with C4orf48 cDNA clone (Cat# [RC227004]) using MegaTran 2.0 (Cat# [TT210002]).