

Product datasheet for **TP327007M**

14-3-3 zeta (YWHAZ) (NM_001135701) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, zeta polypeptide (YWHAZ), transcript variant 5, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227007 representing NM_001135701 Red =Cloning site Green =Tags(s)
	 MDKNELVQKAKLAEQAERYDDMAACMKSVTEQGAELSNEERNLLSVAYKNVVGARRSSWRVSSIEQKTE GAEKKQQMAREYREKIETELRDICNDVLSLLEKFLIPNASQAESKVFLKMKGDYYRYLAEVAAGDDKKG IVDQSQQAYQEAFEISKEMQPHTPIRLGLALNFSVFYIEILNSPEKACSLAKTAFDEAIAELDTLSEES YKDSTLIMQLLRDNLTLWTSDTQGDEAEAGEGGEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	27.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129173
Locus ID:	7534



[View online »](#)

UniProt ID:	P63104 , D0PNI1
Cytogenetics:	8q22.3
RefSeq ORF:	735
Synonyms:	14-3-3-zeta; HEL-S-3; HEL-S-93; HEL4; KCIP-1; POPCHAS; YWHAD
Summary:	This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene. [provided by RefSeq, Oct 2008]
Protein Pathways:	Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis, Pathogenic Escherichia coli infection

Product images:



Coomassie blue staining of purified YWHAZ protein (Cat# [TP327007]). The protein was produced from HEK293T cells transfected with YWHAZ cDNA clone (Cat# [RC227007]) using MegaTran 2.0 (Cat# [TT210002]).