

Product datasheet for **TP327035M**

MYL12B (NM_001144946) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens myosin, light chain 12B, regulatory (MYL12B), transcript variant 4, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC227035 representing NM_001144946
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MSSMFAMFDQSQIQEFKEAFNMIDQNRDGFIDKEDLHDMLASLGKNPTDAYLDAMMNEAPGPINFTMFLT
MFGEKLNKTDPEDEVIRNAFACFDEEATGTIQEDYLRELLTMTMGDRFTDEEVDELYREAPIDKKGNFNIE
FTRILKHGAKDKDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	17.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138418
Locus ID:	103910
UniProt ID:	O14950



[View online »](#)

Cytogenetics:	18p11.31
RefSeq ORF:	462
Synonyms:	MLC-B; MRLC2
Summary:	The activity of nonmuscle myosin II (see MYH9; MIM 160775) is regulated by phosphorylation of a regulatory light chain, such as MRLC2. This phosphorylation results in higher MgATPase activity and the assembly of myosin II filaments (Iwasaki et al., 2001 [PubMed 11942626]). [supplied by OMIM, Mar 2008]
Protein Pathways:	Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton, Tight junction

Product images:



Coomassie blue staining of purified MYL12B protein (Cat# [TP327035]). The protein was produced from HEK293T cells transfected with MYL12B cDNA clone (Cat# [RC227035]) using MegaTran 2.0 (Cat# [TT210002]).