

## Product datasheet for **TP327049M**

### 14-3-3 zeta (YWHAZ) (NM\_001135702) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, zeta polypeptide (YWHAZ), transcript variant 6, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227049 representing NM_001135702 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MDKNELVQKAKLAEQAERYDDMAACMKSVTEQGAELSNEERNLLSVAYKNVVGARRSSWRVSSIEQKTE GAEKKQQMAREYREKIETELRDICNDVLSLLEKFLIPNASQAESKVFLKMKGDYYRYLAEVAAGDDKKG IVDQSQQAYQEAFEISKEMQPHTPIRLGLALNFSVFYEILNSPEKACSLAKTAFDEAIAELDTLSEES YKDSTLIMQLLRDNLTLWTSDTQGDEAEAGEGGEN
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	27.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001129174</a>
Locus ID:	7534



[View online »](#)

UniProt ID: [P63104](#), [D0PNI1](#)

Cytogenetics: 8q22.3


RefSeq ORF: 735

Synonyms: 14-3-3-zeta; HEL-S-3; HEL-S-93; HEL4; KCIP-1; POPCHAS; YWHAD

**Summary:** This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene. [provided by RefSeq, Oct 2008]

**Protein Pathways:** Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis, Pathogenic Escherichia coli infection

### Product images:



Coomassie blue staining of purified YWHAZ protein (Cat# [TP327049]). The protein was produced from HEK293T cells transfected with YWHAZ cDNA clone (Cat# [RC227049]) using MegaTran 2.0 (Cat# [TT210002]).