

Product datasheet for **TP327072L**

SCRN1 (NM_001145513) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human secernin 1 (SCRN1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227072 protein sequence Red =Cloning site Green =Tags(s)

MAAAPPSYCFVAFPPRAKDGLVVF GKNSARPRDEVQEVVYFSAADHEPESKVECTYISIDQVPRTYAIMI
SRPAWLWGAEMGANEHGVCIANEAINTREPAAEIEALLGMDLVRLGLERGETAKEALDVIVSLLEEHGQG
GNYFEDANSCHSFQ SAYLIVDRDEAWVLETIGKYWAAEKVTEGVRCICSQLSLTTKMDAEHPELRSYAQS
QGWWTGEGEFNFSEVFPVEDHLDCGAGKDSLEKQESITVQ TMMNLTLDKASGVCIDSEFFLTTASGVS
VLPQNRSSPCIH YFTGTHDPSRSIFKPFIFVDDVKLVPKTQSPCFGDDDPAKKEPRFQEKPDRRHELYKA
HEWARAIIESDQEQGRKLRSTMLELEKQGLEAMEEILTSSEPLDPAEVGDLFYDCVDTEIKFFK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	46.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001138985
Locus ID:	9805



[View online »](#)

UniProt ID: [Q12765](#), [A0A090N7T9](#)

RefSeq Size: 5286


Cytogenetics: 7p14.3

RefSeq ORF: 1242

Synonyms: SES1

Summary: This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]

Product images:



Coomassie blue staining of purified SCR1 protein (Cat# [TP327072]). The protein was produced from HEK293T cells transfected with SCR1 cDNA clone (Cat# [RC227072]) using MegaTran 2.0 (Cat# [TT210002]).