

Product datasheet for **TP327114L**

EDC3 (NM_001142444) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human enhancer of mRNA decapping 3 homolog (<i>S. cerevisiae</i>) (EDC3), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC227114 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MATDWLGSIVSINCGDSLGVYQGRVSAVDQVSQTISLTRPFHNGVKCLVPEVTFRAGDITELKILEIPGP
GDNQHFGDLHQTELGPSTGAGCQVGINQNGTGKFKKPASSSSAPQNIKRTDVKSDVAVSPQQQCSKS
YVDRHMESLSQSKSFRRRHNSWSSSRHPNQATPKKSLKNGQMKNKDDECFGDDIEEIPDTDFDFEGNL
ALFDKAAVFEEIDTYERRSGTRSRGIPNERPTRYRH DENILESEPIVYRRIIVPHNVSKEFCTDSGLVVP
SISYELHKKLLSVAEKHGLTLERRLEMTGVCASQMALTLLGGPNRLNPKNVHQRPTVALLCGPHVKAQG
ISCGRHLANHDVQVILFLPNFVKMLESITNELSLFSKTQGGQVSSLDLPTSPVDLVINCLDCPENVFLR
DQPWYKAAVAWANQNRAPVLSIDPPVHEVEQGIDAKWSLALGLPLPLGEHAGRIYLCDIGIPQQVFQEVG
INYHSPFGCKFVIPLHSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	55.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001135916
Locus ID:	80153
UniProt ID:	Q96F86
RefSeq Size:	3871
Cytogenetics:	15q24.1
RefSeq ORF:	1524
Synonyms:	hYjeF_N2-15q23; LSM16; MRT50; YJDC; YJEFN2
Summary:	This gene encodes a protein that is important in mRNA degradation. The encoded protein is a component of a decapping complex that promotes efficient removal of the monomethylguanosine (m7G) cap from mRNAs, as part of the 5' to 3' mRNA decay pathway. Mutations in this gene have been identified in human patients with an autosomal recessive form of intellectual disability. [provided by RefSeq, May 2017]
Protein Pathways:	RNA degradation

Product images:



Coomassie blue staining of purified EDC3 protein (Cat# [TP327114]). The protein was produced from HEK293T cells transfected with EDC3 cDNA clone (Cat# [RC227114]) using MegaTran 2.0 (Cat# [TT210002]).