

## Product datasheet for TP327128L

### CCDC124 (NM\_001136203) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 124 (CCDC124), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227128 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MPKKFQGENTKSAAARARRAEAKAAADAKKQKELEDAYWKDDDKHVMRKEQRKEEKEKRRRLDQLERKKET  
QRLLEEEDSKLGGKAPRVATSSKVTRAQIEDTLRRDHQLREAPDTAEKAKSHLEVPLEENVNRRVLEEG  
SVEARTIEDAIAVLSVAEEAADRHPERRMRAAFTAFEEAQLPRLKQENPNMRLSQLKQLLKEWLRSPDN  
PMNQRAVPFNAPK

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	25.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001129675</a>
Locus ID:	115098



[View online »](#)

UniProt ID: [Q96CT7](#), [A0A024R7M8](#)  
RefSeq Size: 1004  
Cytogenetics: 19p13.11  
RefSeq ORF: 669  
Summary: Required for proper progression of late cytokinetic stages.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified CCDC124 protein (Cat# [TP327128]). The protein was produced from HEK293T cells transfected with CCDC124 cDNA clone (Cat# [RC227128]) using MegaTran 2.0 (Cat# [TT210002]).