

Product datasheet for TP327148L

C12orf73 (NM_001135570) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 12 open reading frame 73 (C12orf73), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227148 representing NM_001135570 <div> <div>Red</div>=Cloning site <div>Green</div>=Tags(s) </div> MPAGVPMSTYLMFAASLLAMCAGAEVHRYRDLTIPEIPKRGELKTELLGLKERKHKPQVSQQEEL K <div> <div>TR</div> <div>TRPLEQKLISEEDLAANDILDYKDDDDKV</div> </div>
Tag:	C-Myc/DDK
Predicted MW:	7.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129042
Locus ID:	728568
UniProt ID:	Q69YU5
Cytogenetics:	12q23.3
RefSeq ORF:	213


[View online »](#)

Synonyms: BR; BRAWNIN

Summary: Essential for mitochondrial respiratory chain complex III (CIII) assembly and stability.
[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified C12orf73 protein (Cat# [TP327148]). The protein was produced from HEK293T cells transfected with C12orf73 cDNA clone (Cat# [RC227148]) using MegaTran 2.0 (Cat# [TT210002]).