

Product datasheet for **TP327281L**

Periostin (POSTN) (NM_001135934) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens periostin, osteoblast specific factor (POSTN), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227281 representing NM_001135934 Red =Cloning site Green =Tags(s)

MIPFLPMFSLLLLLLIVNPINANNHYDKILAHSRIRGRDQGPNCALQQILGTTKYYFSTCKNWWYKKSICG
QKTTVLYECCPGYMRMEGMKGCPAVLPIDHVVYGLGIVGATTTQRYSDASKLREEIEGKGSFTYFAPSNE
AWDNLDSDIRRGLESNVNVELLNALHSHMINKRMLTKDLKNGMIIPSMYNNLGLFINHYPNGVVTVNCAR
IIHGNQIATNGVHVHIDRVLTQIGTSIQDFIEAEDDLSSFRAAAITSDILEALGRDGHFTLFAPTNEAFE
KLPRGVLERIMGDKVASEALMKYHILNTLQCSSEIMGGAVFETLEGNTIEIGCDGDSITVNGIKMVNKKD
IVTNNGVIIHLIDQVLIPDSAKQVIELAGKQQTTFDLVAQLGLASALRPDGEYLLAPVNNAFSDDTLSM
DQRLKLLQNHILKVKVGLNELYNGQILETIGGKQLRVFVYRTAVCIENSCMEKGSKQGRNGAIHFIRE
IIPAEKSLHEKQDKRFSTFLSLEAADLKELLTQPGDWTLFVPTNDAFKGMTSEEKILIRDKNALQ
NIILYHLTPGVFIGKGFEPGVTNLIKTTQGSKIFLKEVNDLLVNELKSKESDIMTTNGVIHVVDKLLYP
ADTPVGNQDQLEILNKLIKYIQIKFVRGSTFKEIPVTYKPIIKKYTKIIDGVPVEITEKETREERIITG
PEIKYTRISTGGGETEETLKKLLQEEVTKVTKFIEGGDGHLEFEEIKRLLQGDTPVRKLQANKKVQGSR
RRLREGRSQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	86.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001129406](#)

Locus ID: 10631

UniProt ID: [Q15063](#), [A0A024RDT5](#)

Cytogenetics: 13q13.3


RefSeq ORF: 2337

Synonyms: OSF-2; OSF2; PDLPOSTN; PN

Summary: This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified POSTN protein (Cat# [TP327281]). The protein was produced from HEK293T cells transfected with POSTN cDNA clone (Cat# [RC227281]) using MegaTran 2.0 (Cat# [TT210002]).