

Product datasheet for **TP327305M**

SUPT20HL2 (NM_001136233) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens family with sequence similarity 48, member B2 (FAM48B2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227305 representing NM_001136233 Red =Cloning site Green =Tags(s)

MDRDLEQALDRTENITEIAQRRPRRRYSRAGKTLQEKLVDIYVEECGKEPEDPQELRSNVNLEKLVR
RESLPCLLVNLYPGNQGYSVMLQREDGSFAETIRLPYEERALLDYLDAAELPPALGDVLDKASVNIHSG
CVIVEVRDYRQSSNMQPPGYQSRHILLRPTMQTLADEVKTMTRDGEKWSQEDKFPLESQILATAEPLCL
DPSVAVACTANRLLYNKQKMNTDPMEQCLQRYSWPSVKPQQEQSDCPPPELRVSTSGQKEERKVGQP
CE
LNITKAGSCVDTWKGRPCDLAVPSEVDVEKLAKGYQSVTAADPQLPVWPAQEVEDPFRHAWAEAGCQAW
DT
KPNIMQSFNDPLLCKGIRPRKKARQKSQKSPWQPFDDHSACLRPGSETDAGRAVSQAQESVQSKVKG
G
KMSHSSSGPASVSQLSSWKTPEQPDVWVQSSVSGKGEKHPPPRRTQLPSSSGKISSGNSFPPQQAGSPL
K
RPFPAAPAVAAAAPAPAPAPAAAPALAAAATAAAAGGAAPSHSQKPSVPLIKASRRRPAAGRPTRFVKI
APAIQVRTGSTGLKATNVEGPVRGAQVLGCSFKPVQAPGSGAPAPAGISGSLQSSGGPLPDARPGAVQA
SSPAPLQFFLNTPEGLRPLTLQVPQGWAVLTGPQQQSHQLVSLQLLQQTAAHPPQPGPQGSTLGLST
QG
QAFPAQQLLNVNLTGAGSGLQPQPQAAVLSLLGSAQVPQQGVQLPFVLGQQPQPLLLLQPQPQPQQI
QLQ
TQPLRVLQQPVFLATGAVQIVQPHPGVQAGSQLVGQRKGGKPTPPAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	87.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method



[View online »](#)

Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001129705</u>
Locus ID:	170067
UniProt ID:	<u>P0C7V6</u>
Cytogenetics:	Xp22.11
RefSeq ORF:	2451
Synonyms:	FAM48B1; FAM48B2

Product images:



Coomassie blue staining of purified SUPT20HL2 protein (Cat# [TP327305]). The protein was produced from HEK293T cells transfected with SUPT20HL2 cDNA clone (Cat# [RC227305]) using MegaTran 2.0 (Cat# [TT210002]).