

Product datasheet for TP327317M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

FGF1 (NM_001144892) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens fibroblast growth factor 1 (acidic) (FGF1),

transcript variant 4, 100 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC227317 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MAEGEITTFTALTEKFNLPPGNYKKPKLLYCSNGGHFLRILPDGTVDGTRDRSDQHIQLQLSAESVGEVY IKSTETGQYLAMDTDGLLYGSQTPNEECLFLERLEENHYNTYISKKHAEKNWFVGLKKNGSCKRGPRTHY

GQKAILFLPLPVSSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 17.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001138364

Locus ID: 2246

UniProt ID: P05230





RefSeq Size: 3682

Cytogenetics: 5q31.3 RefSeq ORF: 465

Synonyms: AFGF; ECGF; ECGF-beta; ECGFA; ECGFB; FGF-1; FGF-alpha; FGFA; GLIO703; HBGF-1; HBGF1

Summary: The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF

family members possess broad mitogenic and cell survival activities, and are involved in a

variety of biological processes, including embryonic development, cell growth,


morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis. Multiple alternatively spliced variants encoding different

isoforms have been described. [provided by RefSeq, Jan 2009]

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton

Product images:



Coomassie blue staining of purified FGF1 protein (Cat# [TP327317]). The protein was produced from HEK293T cells transfected with FGF1 cDNA clone (Cat# [RC227317]) using MegaTran 2.0 (Cat# [TT210002]).