

## Product datasheet for **TP327319L**

### TTC39A (NM\_001144832) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tetratricopeptide repeat domain 39A (TTC39A), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227319 representing NM_001144832 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAFLTTWKEAASGKSDRRTPESSLHEALDQCMTALDLFLTNQFSEALSYLKPRTKESMYHSLTYATILEM  
QAMMTFDPQDILLAGNMMKEAQMCLCQRHRRKSSVTDSFSSLVNRPTLGQFTEEEIHAEVCYAECLLQRAA  
LTFLQDENMVSFIKGGIKVRNSYQTYKELDSLQSSQYCKGENHPHFEGGVKLGVGAFNLTLMSLPTRIL  
RLLEFVGFSGNKDYGLLQLEEGASGHSFRSVLCVMLLLCYHTFLTFVLGTGNVNIEEAELLKPYLNRYP  
KGAIFLFFAGRIEVIKGNIDAAIRRFEECCEAQQHWKQFHHCYWELMWCFYKQWKMSYFYADLLSKE  
NCWSKATYIYMKAAYLSMFGKEDHKPFGDDEVELFRAVPGLKLLKIAGKSLPTEKFAIRKSRRYFSSNPIS  
LPVPALEMMYIWNQYAVIGKQPKLTDGILEITKAEEMLEKGPENEYSVDDECLVKLLKGLCLKYLGRVQ  
EAEENFRSISANEKIKYDHYLIPNALLELALLLMEQDRNEEAIKLLESKQNYKNYSMESRTHFRIQAA  
TLQAKSSLENSRSMVSSVSL


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	66.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001138304</a>
<b>Locus ID:</b>	22996
<b>UniProt ID:</b>	<a href="#">Q5SRH9</a>
<b>Cytogenetics:</b>	1p32.3
<b>RefSeq ORF:</b>	1743
<b>Synonyms:</b>	C1orf34; DEME-6

**Product images:**

Coomassie blue staining of purified TTC39A protein (Cat# [TP327319]). The protein was produced from HEK293T cells transfected with TTC39A cDNA clone (Cat# [RC227319]) using MegaTran 2.0 (Cat# [TT210002]).