

Product datasheet for TP327328

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Annexin A2 (ANXA2) (NM_001136015) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human annexin A2 (ANXA2), transcript variant 4, 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC227328 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MSTVHEILCKLSLEGDHSTPPSAYGSVKAYTNFDAERDALNIETAIKTKGVDEVTIVNILTNRSNAQRQD IAFAYQRRTKKELASALKSALSGHLETVILGLLKTPAQYDASELKASMKGLGTDEDSLIEIICSRTNQEL QEINRVYKEMYKTDLEKDIISDTSGDFRKLMVALAKGRRAEDGSVIDYELIDQDARDLYDAGVKRKGTDV PKWISIMTERSVPHLQKVFDRYKSYSPYDMLESIRKEVKGDLENAFLNLVQCIQNKPLYFADRLYDSMKG

KGTRDKVLIRIMVSRSEVDMLKIRSEFKRKYGKSLYYYIQQDTKGDYQKALLYLCGGDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 38.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001129487

Locus ID: 302



UniProt ID: <u>P07355</u>, <u>A0A024R5Z7</u>

RefSeq Size: 1665 Cytogenetics: 15q22.2 RefSeq ORF: 1017

Synonyms: ANX2; ANX2L4; CAL1H; HEL-S-270; LIP2; LPC2D; P36; PAP-IV


Summary: This gene encodes a member of the annexin family. Members of this calcium-dependent

phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. Annexin A2 expression has been found to correlate with resistance to treatment against various cancer forms. [provided by

RefSeq, Dec 2019]

Protein Families: Druggable Genome, Secreted Protein, Stem cell - Pluripotency

Product images:



Coomassie blue staining of purified ANXA2 protein (Cat# TP327328). The protein was produced from HEK293T cells transfected with ANXA2 cDNA clone (Cat# [RC227328]) using MegaTran 2.0 (Cat# [TT210002]).