

Product datasheet for **TP327429**

RFPL4A (NM_001145014) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ret finger protein-like 4A (RFPL4A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227429 representing NM_001145014 Red =Cloning site Green =Tags(s) MAEHFKQIIRCPVCLKDLEEAVQLKCGYACCLQCLNSLQKEPDGEGLLCRFCSVVSQKDDIKPKYKLRL VSIKELEPKLKSVLTMNPRMRKFQVDMTFDVTANNYLIISEDLSFRSGDLSQNRKEQAERFDALCV LGTPRFTSGRHYWEVDVGTSQVWDVGCKESVNRQGKIVLSSEHGFLTVGCREGKVFAASTVPMTPLVV S PQLHRVGIFLDVGMRSIAFYNVSDGCHIYTFIEIPVCEPWRPFFAHKRGSQDDQSILSICSVINPSAASA PVSSEGK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001138486</u>
Locus ID:	342931



[View online »](#)

UniProt ID: [A6NLU0](#)
 Cytogenetics: 19q13.42
 RefSeq ORF: 861
 Synonyms: RFPL4; RNF210
 Protein Families: Druggable Genome

Product images:


Coomassie blue staining of purified RFPL4A protein (Cat# TP327429). The protein was produced from HEK293T cells transfected with RFPL4A cDNA clone (Cat# [RC227429]) using MegaTran 2.0 (Cat# [TT210002]).