

## Product datasheet for **TP327432**

### TTC39C (NM\_001135993) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins  
**Description:** Recombinant protein of human tetratricopeptide repeat domain 39C (TTC39C), transcript variant 1, 20 µg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA Clone or AA Sequence:** >RC227432 representing NM\_001135993  
**Red**=Cloning site **Green**=Tags(s)

MAGSEQRRPRRRDDGSDSAAAAAAPLQDAELALAGINMLLNNGFRES DQLFKQYRNHSPLMSFGASFVS  
FLNAMMTFEEEEKMLACDDLKTEKLCSEEEAGVIETIKNKIKKNVDVRKSAPSMVDRLQRQIIIADCQV  
YLAVLSFVKQELSAYIKGGWILRKAWKIYNKCYLDINALQELYQKKLTEESLTS DAANDNHIVAEGVSEE  
SLNRLKGAVSFGYGLFHLCISMVPPNLLKIINLLGFPGDRLQGLSSLMYASESKDMKAPLATLALLWYHT  
VWRPFFALDGSDNKAGLDEAKEILLKKEAAYPNSSLFMFFKGRIQRLEQCINSALTSFHTALELAVDQRE  
IQHVCLYEIGWCSMIELNFKDAFDSFERLKNESRWSQCYAYLTAVCQGATGDVDGAQIVFKEVQKLFKR  
KNNQIEQFSVKKAERFRKQTPTKALCVLASIEVLYLWKALPNCSFPNLQRMSQACHEVDDSSVGLKYLL  
LGAIHKCLGNSEDAVQYFQRAVKDELRCRQNNLYVQPYACYELGCLLLDKPETVGRGRALLQAKEDFSGY  
DFENRLHVRIHAALASLRELVPQ


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK  
**Predicted MW:** 65.7 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  
**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.  
**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.  
**Storage:** Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001129465</a>
<b>Locus ID:</b>	125488
<b>UniProt ID:</b>	<a href="#">Q8N584</a>
<b>Cytogenetics:</b>	18q11.2
<b>RefSeq ORF:</b>	1749
<b>Synonyms:</b>	C18orf17; HsT2697

**Product images:**

Coomassie blue staining of purified TTC39C protein (Cat# TP327432). The protein was produced from HEK293T cells transfected with TTC39C cDNA clone (Cat# [RC227432]) using MegaTran 2.0 (Cat# [TT210002]).