

Product datasheet for **TP327432M**

TTC39C (NM_001135993) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tetratricopeptide repeat domain 39C (TTC39C), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227432 representing NM_001135993 Red =Cloning site Green =Tags(s)

MAGSEQQRPRRRDDGDSDAAAAAAALPQDAELALAGINMLLNNGFRES DQLFKQYRNHSPLMSFGASFVS
FLNAMMTFEEEEKMLACDDLKTEKLCSEEEAGVIETIKNKIKKNVDVRKSAPSMVDRLQRQIIIADCQV
YLAVLSFVKQELSAYIKGGWILRKAWKIYNKCYLDINALQELYQKKLTEESLTSDAANDNHIVAEGVSEE
SLNRLKGAVSFGYGLFHLCSMVPPNLLKIINLLGFPGDRLQGLSSLMYASESKDMKAPLATLALLWYHT
VWRPFFALDGSDNKAGLDEAKEILLKKEAAYPNSSLFMFFKGRIQRLECQINSALTSFHTALELAVDQRE
IQHVCLYEIGWCSMIELNFKDAFDSFERLKNESRWSQCYAYLTAVCQGATGDVDGAQIVFKEVQKLFKR
KNNQIEQFSVKKAERFRKQTPTKALCVLASIEVLYLWKALPNCSFPNLQRMSQACHEVDDSSVGLKYLL
LGAHKCLGNSEDAVQYFQRAVKDELRCRQNNLYVQPYACYELGCLLLDKPETVGRGRALLLQAKEDFSGY
DFENRLHVRIHAALASLRELVPQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	65.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001129465
Locus ID:	125488
UniProt ID:	Q8N584
Cytogenetics:	18q11.2
RefSeq ORF:	1749
Synonyms:	C18orf17; HsT2697

Product images:

Coomassie blue staining of purified TTC39C protein (Cat# [TP327432]). The protein was produced from HEK293T cells transfected with TTC39C cDNA clone (Cat# [RC227432]) using MegaTran 2.0 (Cat# [TT210002]).