

Product datasheet for **TP327595**

ECSIT (NM_001142465) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ECSIT homolog (Drosophila) (ECSIT), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC227595 representing NM_001142465 Red =Cloning site Green =Tags(s) MSWVQATLLARGLCRAWGGTCGAALTGTSISQVPLPKDSTGAADPPQPHIVGIQSPDQQAALARHNPARP VFVEGPFSLWLRNKCVYYHILRADLLPPEEREVEETPEEWNLYYPMQLDLEYVRSGWDNYEFDINEVEEG PVFAMCMAGAHDQATMAKWIQGLQETNPTLAQIPVVFRLAGSTRELQTSAGLEEPPLPEDHQEEDDLNQL RQQGGQS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001135937
Locus ID:	51295
UniProt ID:	Q9BQ95



[View online »](#)

Cytogenetics: 19p13.2


RefSeq ORF: 651

Synonyms: SITPEC

Summary: Adapter protein of the Toll-like and IL-1 receptor signaling pathway that is involved in the activation of NF-kappa-B via MAP3K1. Promotes proteolytic activation of MAP3K1. Involved in the BMP signaling pathway. Required for normal embryonic development (By similarity). [UniProtKB/Swiss-Prot Function]

Protein Pathways: MAPK signaling pathway

Product images:



Coomassie blue staining of purified ECSIT protein (Cat# TP327595). The protein was produced from HEK293T cells transfected with ECSIT cDNA clone (Cat# [RC227595]) using MegaTran 2.0 (Cat# [TT210002]).